

ш

ш

工

S

◁

<

4

0

Z

工

0

ш

## TECHNICAL DATA SHEET

## TYPE 6 STEAM CHEMICAL INDICATOR (3,5 min)

PRODUCT CODE 4AC6

4A Medical Steam Indicator is clasified as Type 6 according to ISO 11140-1 **RELATED STANDARDS** 

4A Medical Steam Indicators are cycle verification indicators that confirm the INTENDED USE

presence or absence of specific time and temperature parameters during a cycle. The products allows the end user to identify incorrect packaging, load problems and

sterilization failures.

4A Medical Steam Indicators are designed operate at steam cycle 134°C - 3,5 **OPERATION CONDITIONS:** 

minutes or 121 °C - 15 minutes operation conditions.

CHEMICAL COMPOSITION: The sterilization indicator ink used for printing contains following components in ppm levels.

> Water 35% (CAS # 7732-18-5) (CAS # 123-42-2) Diacetone Alcohol 3% **Proprietary Mixture** 62% (CAS # N/A)

> > **Before Sterilization After Sterilization**

**PRODUCT** 

Each box contains 250 pieces indicator printed carton sheets. The indicator is temperature, time and steam concentration sensitive. **SPECIFICATIONS** 

STORAGE AND DISPOSAL The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and

defective products are not used. Shelf life is 5 years after prodcution date. After

sterilization indicator keep sterilization colour for 6 month.

100 boxes in cardboard. **PACKAGING**