

TECHNICAL DATA SHEET

ETHYLEN OXIDE CHEMICAL INDICATOR

TECHNICAL DATA SHEET

-
- PRODUCT CODE** : 4AEOI
- RELATED STANDARDS** : 4A Medical Ethylen Oxide Chemical Indicator is classified as Type 5 according to ISO 11140-1 standard.
- INTENDED USE** : 4A Medical Ethylen Oxide is chemical indicator which gives information about success of sterilization by using within sterilization package together with ethylene oxide.
- OPERATION CONDITIONS** : It is used to evaluate the process of Ethylene Oxide Gas sterilization. 4A indicators which were produced according to ISO 11140 are water based and they do not include toxic substances.

Before Sterilization

After Sterilization



- PRODUCT SPECIFICATIONS** : Each box contains 500 pieces indicator printed carton sheets. The indicator is ethylen oxide concentration sensitive.
- STORAGE AND DISPOSAL** : The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.
- PACKAGING** : 20 packs in cardboard.