

Ш

S

4

CA

I

0

Ш

## **TECHNICAL DATA SHEET**

## **CONTAINER FILTER**

PRODUCT CODE : 4ACF

**RELATED STANDARDS**: 4A Medical Container Filter's Indicator is clasified as Type1 according to

ISO 11140-1 standard.

INTENDED USE : 4A Medical Container filter is used in sterilization containers as a filter and

indicator is to control sterilization process. 4A Medical has round and

rectangular container filter.

**OPERATION CONDITIONS**: 4A Medical Container Filter are designed to operate at steam cycle. Filter is

for single use.

CHEMICAL COMPOSITION: It is made of white crepe paper and non toxic sterilization ink. The

sterilization indicator ink used for printing contains following components

in ppm levels

 Water
 35%
 (CAS # 7732-18-5)

 Diacetone Alcohol
 3%
 (CAS # 123-42-2)

 Proprietary Mixture
 62%
 (CAS # N/A)

Contains no substance specified as hazardous, which should be taken into

account according to EC directives.

Before Sterilization After Sterilization

**PRODUCT** 

SPECIFICATIONS

Each pack contains 100 pieces. The indicator is temperature and steam

concentration sensitive.

STORAGE AND DISPOSAL

The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact

: with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization

colour for 6 month.

PACKAGING : Packed in cardboxes.