

Ш

Ш

工

S

<

4 

4

0

Z

I

0

ш

## **TECHNICAL DATA SHEET**

## INSTRUMENT PROTECTOR

PRODUCT CODE

RELATED STANDARDS 4A Medical Instrument Protector's Indicator is clasified as Type1 according to

ISO 11140-1 standard.

4AIP

4A Medical Instrument Protector keeps the tips of instruments from penetrating INTENDED USE

the sterile pouch. It holds the instrument in an open position for proper

sterilization. It is for use inside a single peel pouch.

**OPERATION CONDITIONS**: 4A Medical Instrument Protectors are designed to operate at steam cycle.

CHEMICAL COMPOSITION: The sterilization indicator ink used for printing contains following components

in ppm levels

Water 35% (CAS # 7732-18-5) Diacetone Alcohol 3% (CAS # 123-42-2) 62% (CAS # N/A) Proprietary Mixture

Contains no substance specified as hazardous, which should be taken into

account according to EC directives.

**Before Sterilization** After Sterilization

**PRODUCT** 

Each pack contains 100 pieces. The indicator is temperature, time and steam SPECIFICATIONS concentration sensitive. Standard dimensions are 50mm x 125mm and 75mm x

150mm. Any other dimension can me made on request.

STORAGE AND DISPOSAL The expired products are not used. It's kept in the places with room temperature

without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date.

After sterilization indicator keep sterilization colour for 6 month.

**PACKAGING** Packed in cardboxes.