













## TECHNICAL DATA SHEET

### PLASMA DOCUMENTATION LABEL

TECHNICAL DATA SHEET

<b>PRODUCT CODE</b>	:	4APDL						
<b>RELATED STANDARDS</b>	:	4A Medical Plasma Documentation Label is classified as Type 1 Indicator according to the ISO 11140-1 standard.						
<b>INTENDED USE</b>	:	4A Medical Plasma Documentation Labels are allow the end user to identify loading problems and sterilization failures.						
<b>OPERATION CONDITIONS</b>	:	4A Medical Plasma Documentation Labels are designed to operate at hydrogen peroxide cycle.  Contains no substance classified as hazardous, which should be taken into account according to EC directives.						
		<table border="0" style="width: 100%;"> <tr> <td style="text-align: center;"><u>Before Sterilization</u></td> <td style="text-align: center;"><u>After Sterilization</u></td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> </table>	<u>Before Sterilization</u>	<u>After Sterilization</u>				
<u>Before Sterilization</u>	<u>After Sterilization</u>							
								
								
<b>PRODUCT SPECIFICATIONS</b>	:	Each roll contains 400 pieces of self adhesive stickers (double adhesive) with indicator. The indicator is printed on the sticker and it is temperature, time and steam concentration sensitive.						
<b>STORAGE AND DISPOSAL</b>	:	The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.						
<b>PACKAGING</b>	:	100 rolls in one cardboard.						