













TECHNICAL DATA SHEET

HYDROGEN PEROXIDE (PLASMA) PLASMA CHEMIAL INDICATOR

TECHNICAL DATA SHEET

PRODUCT CODE	:	4API						
RELATED STANDARDS	:	Product is classified as Type 4 Indicator according to the ISO 11140-1 standard.						
INTENDED USE	:	Plasma sterilization indicator is used in low temperature hydrogen peroxide plasma sterilization process instructions.						
OPERATION CONDITIONS	:	It is used to evaluate the process of HP Plasma sterilization. Correct color change shows whether item exposed to sterilisation process or not. HP Plasma indicator is water-based and non-toxic chemical indicator.						
		<table border="0" style="width: 100%;"> <tr> <td style="text-align: center;"><u>Before Sterilization</u></td> <td style="text-align: center;"><u>After Sterilization</u></td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> </table>	<u>Before Sterilization</u>	<u>After Sterilization</u>				
<u>Before Sterilization</u>	<u>After Sterilization</u>							
								
								
PRODUCT SPECIFICATIONS	:	Each box contains 500 pieces indicator printed carton sheets. The indicator is hydrogen peroxide concentration sensitive.						
STORAGE AND DISPOSAL	:	The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.						
PACKAGING	:	20 packs in cardboard.						