

TECHNICAL DATA SHEET

PLS WASHER TEST INDICATOR

PRODUCT CODE	:	4APLS
RELATED STANDARDS	:	4A Medical PLS Washer Test is under the scope of ISO 13485 and 9001 standards.
INTENDED USE	:	PLS Washer test is used to control whether cleaning process is successful or not.
OPERATION CONDITIONS	:	Surgical devices are placed into the wire basket before washing process. PLS Washer test has stainless steel frame and it is placed in it. PLS Washer test is placed into the wire basket. It is sent to the washing process. If there is any residue on the washer test after washing process, it means that washing process is not successful and it should be repeated. If there is no residue it means washing process is successful.
PRODUCT SPECIFICATIONS	:	Each chemical washer test is disposable. (single use only). Because of its design it is easy for documentation.
STORAGE AND DISPOSAL	:	Do not use damaged products. Keep products at room temperature 25°C.
PACKAGING	:	100 pcs in a box. 5000 pcs in a cardboard.