



TECHNICAL DATA SHEET

SONI CONTROL

PRODUCT CODE : 4ASC

RELATED STANDARDS : 4A Medical Soni Control is under the scope of ISO 13485 and 9001

standards.

INTENDED USE : Sonicontrol product is used for testing sonication performance of ultrasonic

washers. This product contains pH sensitive indicator materials which change its color according to changing pH conditions during sonication

process.

OPERATION CONDITIONS: The amounts of indicator substances were carefully adjusted at desired level

to make color change significant during the process. Depending on the temperature and environmental conditions, pH of the solution can shift a little bit and therefore small color changes (oscillations) between green and blue starting colors will be able to occur in solution at these amount levels for indicator materials. However this situation cannot affect performance of

the product.

STORAGE AND DISPOSAL Store at temperatures between 2-25°C. Do not freeze. Keep away from

light, heat and chemicals.

PACKAGING : 24 pcs in a box.