













## TECHNICAL DATA SHEET

### STEAM DOCUMENTATION LABEL

TECHNICAL DATA SHEET

<b>PRODUCT CODE</b>	:	4ASDL															
<b>RELATED STANDARDS</b>	:	4A Medical Steam Documentation Label is classified as Type 1 Indicator according to the ISO 11140-1 standard.															
<b>INTENDED USE</b>	:	4A Medical Steam Documentation Labels are allow the end user to identify loading problems and sterilization failures.															
<b>OPERATION CONDITIONS</b>	:	4A Medical Steam Documentation Labels are designed to operate at steam cycle.															
<b>CHEMICAL COMPOSITION</b>	:	<p>The sterilization indicator ink used for printing contains following components in ppm levels</p> <table border="0" style="margin-left: 40px;"> <tr> <td>Water</td> <td style="text-align: right;">35%</td> <td style="text-align: right;">(CAS # 7732-18-5)</td> </tr> <tr> <td>Diacetone Alcohol</td> <td style="text-align: right;">3%</td> <td style="text-align: right;">(CAS # 123-42-2)</td> </tr> <tr> <td>Proprietary Mixture</td> <td style="text-align: right;">62%</td> <td style="text-align: right;">(CAS # N/A)</td> </tr> </table> <p>Contains no substance classified as hazardous, which should be taken into account according to EC directives.</p> <table border="0" style="margin-left: 40px; width: 100%;"> <thead> <tr> <th style="text-align: center;"><u>Before Sterilization</u></th> <th style="text-align: center;"><u>After Sterilization</u></th> </tr> </thead> <tbody> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> <tr> <td style="text-align: center;"></td> <td style="text-align: center;"></td> </tr> </tbody> </table>	Water	35%	(CAS # 7732-18-5)	Diacetone Alcohol	3%	(CAS # 123-42-2)	Proprietary Mixture	62%	(CAS # N/A)	<u>Before Sterilization</u>	<u>After Sterilization</u>				
Water	35%	(CAS # 7732-18-5)															
Diacetone Alcohol	3%	(CAS # 123-42-2)															
Proprietary Mixture	62%	(CAS # N/A)															
<u>Before Sterilization</u>	<u>After Sterilization</u>																
																	
																	
<b>PRODUCT SPECIFICATIONS</b>	:	Each roll contains 500 pieces of self adhesive stickers (double adhesive) with indicator. The indicator is printed on the sticker and it is temperature, time and steam concentration sensitive.															
<b>STORAGE AND DISPOSAL</b>	:	The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.															
<b>PACKAGING</b>	:	150 rolls in one cardboard.															